**Comfy**Splints® We Make The World Comfy

# Pediatrics

### **Comfortable** Effective Splints

for even the smallest users

(800) 582-5332 www.comfysplints.com



### **Benefits & Features of Comfy Splints**

- Soft and cuddly
- Multilayered for comfort •
- Easy to use •
- Patented "bend-to-fit" customizable frame •
- Reinforced spine can be adjusted and readjusted without tools .
- Removable, machine-washable covers •
- Static, progressive and dynamic models available •
- Latex-free, antimicrobial and moisture-wicking materials for skin breakdown prevention



## Indications Comfy Splints<sup>®</sup>

- Sprains / Strains / Fractures
- Neurological impairments
- Contractures
- Spasticity
- CP

# Loved by Doctors, Therapists, patients and Families **Benefits & Features** of Comfyprene Splints<sup>®</sup>

- Patented bend-to-fit customizable frame
- Soft neoprene (without removable covers)
- Offers the advantages of our Comfy line, in a more slim, discrete form
- Easily hand washable with soap and water
- Comes with hygiene liners for extra protection
- Straps can be positioned anywhere for optimal fit
- Four color options (Light Blue, Purple, Navy and Camouflage)
- Ideal for sports injuries, cast removal, sports medicine and more



Comfyprene color options shown: Light Blue, Purple, Navy and Camouflage Available in all sizes

### **Comfy Hand Orthoses**

- Support and positioning of the weak or deformed hands
- Does not stretch thenar eminence
- Resists strong flexion synergy
- Super soft straps
- Kid-friendly color options
- Static, progressive and dynamic models available
- Arm, wrist and fingers are independently adjustable
- HCPCS (PDAC) billing codes available



### Comfy Elbow Orthoses

- Adjustable and ambidextrous cuffs
- Positions the flaccid elbow and arm in functional position
- Provides gradual extension of the elbow
- Bifurcated center straps
- Helps restrict elbow flexion pattern problems and deformities or neuromuscular impairment
- Super comfortable to wear
- Static, progressive and dynamic models available
- HCPCS (PDAC) billing codes available

Did you know that any joint (knee, elbow or wrist) can be turned into a dynamic moveable one by adding our patented, simple to use goniometer mechanism?

> Comfyprene Elbow Orthosis in Navy

Comfy Knee Orthosis in Pediatric Large

Comfy Elbow Orthosis in Pediatric Small



Goniometers allow for a set range of motion (ROM) or desired static set position in 5° increments. Spring Loaded Goniometer option has additional spring load that gently pushes extremity to selected degree of extension. (Can be special ordered to increase flexion.)

> Comfyprene Knee Orthosis in Light Blue color

Comfy Pediatric Boot

### Comfy Ankle & Foot Orthoses

- Universal sizing fits most
- Adjustable toe post
- Anti hip-rotation bar
- Ambulating/transfer pad available
- Static, progressive and dynamic models available
- HCPCS (PDAC) billing codes available

### Comfy Knee Orthoses

- Adjustable and ambidextrous cuffs
- Soft edges prevent pressure areas
- Bifurcated center straps
- Stabilize/immobilize joints in desired position
- Gradually stretch to increase range of motion
- Provide strength for flaccid or weak extremity
- Static, progressive and dynamic models available
- HCPCS (PDAC) billing codes available

Comfy Pediatric Hands (PH101; PHT101; POPH101; PDH101; PGH101; PF101)							
Size	Size Wrist to Tip of Middle Finger Total Splint Length Average Age						
Pediatric Large	4 ½" – 6"	9 1⁄2"	9 and up				
Pediatric Medium 3" - 4 ½" 8" 4-9 yrs. old							
Pediatric Small 2 ½" – 3 ½" 7" 1-4 years old							

Comfyprene Hands (H-CP-P; HT-CP-P; OPH-CP-P; GH-CP-P; F-CP-P)						
Size	Size Wrist to Tip of Middle Finger Total Length Average Age					
Pediatric Large	4 ½" – 6"	9"	9 and up			
Pediatric Medium 3" - 4 ½" 7" 4-9 yrs. old						
Pediatric Small 2 ½" – 3 ½" 6 ½" 1-4 years old						

Comfy Adjustable Cone Hand (PACH101)						
Size Wrist to Tip of Middle Finger Total Length Average Age						
Pediatric Large	4 1/2" - 6"	9 1⁄2"	9 and up			
Pediatric Medium	3" – 4 ½"	8″	4-9 yrs. old			
Pediatric Small	2 1/2" - 3 1/2"	7"	1-4 years old			

Comfyprene Separate Finger Hand (HSF-CP-P)					
Size Wrist to Tip of Middle Finger Total Length Average Age					
Pediatric Large	4 1⁄2″ – 6″	10"	9 and up		

Comfy Elbow-Hand (PEH101)					
Size Bicep Circumference Splint Length Average Age					
Pediatric Large	8"-10"	11"-14"	9 and up		
Pediatric Medium	6"-8"	9"-13"	4-9 yrs. old		
Pediatric Small	4"-6"	7"-12"	1-4 years old		

Comfy Knees (PK101; PGK101; PSGK101) Comfyprene Knees (K-CP-P; GK-CP-P; SGK-CP-P)					
Size Knee Crease Knee Crease towards Ankle towards Thigh Total Splint Length Average Age					
Pediatric Large 5 ¼" 5 ¼" 10 ½" 9 and up					
Pediatric Medium	4 ½"	4 ½"	9″	4-9 yrs. old	
Pediatric Small	3 ¾"	3 ¾"	7 ½″	1-4 years old	

Comfy Elbows (PE101;, PGE101; PSGE101) Comfyprene Elbows ( E-CP-P; GE-CP-P; SGE-CP-P)				
Size Elbow Crease Elbow Crease towards Total Splint Average towards Wrist Shoulder Length Age				
Pediatric Large	4 ½"	4 1⁄2″	9	9 and up
Pediatric Medium	3 ¾	3 3/4	7 ½"	4-9 yrs. old
Pediatric Small	3 ¼"	3 ¼"	6 ½"	1-4 yrs. old

Lenjoy Medical Engineering, Inc. is committed to providing innovative, effective and cutting edge products that are comfortable, easy to use and affordable. Our team of Occupational Therapists, M.D.'s, nurses and engineers are continuously working to meet the needs of those in need of our products. In our continuing efforts, we will make reasonable and necessary adjustments upon your request. We know that you will soon see why therapists, doctors, patients and families have come to love our products.

#### **Returns policy:**

We stand by our products and will exchange any defective products (in the unlikely event).

We ask that you call the Lenjoy offices (800-582-5332) before any return to obtain the necessary return number and repayment authorization. We find that people are happy to talk with our professional staff as correct product ordering, not product efficacy or quality, is often times the main barrier.

- Returns of unused (fully re-saleable) items may be exchanged or refunded within 30 days of invoice
- Returned items will not be accepted and credit will not be issued without prior authorization by Lenjoy Medical Engineering, Inc.
- All returns must have an original packing slip, order number or P.O. number
- All returned items are subject to review before credit is issued

- Custom adaptations, special order items, used or altered items will not be accepted for credit.

- Purchaser is responsible for shipping and handling costs of returned items.
- Authorized returns must be sent freight prepaid.
- Any product returned is subject to a 15% restocking charge.
- No returns can be accepted after 90 days past invoice date.

#### Warranty:

LME warrants it's products for single use for 90 days. Failure to follow device instructions voids warranty.

#### **Prices and Terms:**

Terms are net 30 days (unless otherwise agreed). We accept Mastercard, Visa, Discover Card, Amex and check payments.

#### Shipping:

All shipments are FOB Lenjoy Medical Engineering, Inc., Gardena, California. All Shipping and handling charges are paid by the customer and are added to the final invoice. Shipments will be sent via UPS ground unless 3 day air, 2 day air or overnight delivery is specified. Orders received by 3:00 pm PST will be shipped same day. Custom orders will be shipped ASAP. If there are any delays or backorders, we will quickly notify customer.

#### Collection:

If collection proceedings are necessary in order to obtain payment on an invoice, the customer will be responsible for paying all collection costs and legal fees. All necessary collections will be adjudicated in the State of California.

#### Please note:

Comfy Splints are for single person use only, require a physician's order for reimbursement, and should be administered by personnel trained in the use and care of splints. Comfy Splint care instructions should be followed at all times. All applied splints should be frequently checked and adjusted to prevent areas of pressure or incorrect use.



Featured on cover: Comfyprene Knee and Hand Separate Finger in Light Blue Featured below: Comfyprene Hand Thumb with finger separator in Purple

> Lenjoy Medical Engineering, Inc. 13112 Crenshaw Blvd. Gardena, CA 90249 Phone (800) 582-5332 • Fax (310) 353-2484 www.comfysplints.com