



CLINICALLY SUPERIOR!

Hands / Wrists / Fingers

* Right/Left ** Ambidextrous

Sample Indication	Recommended Orthosis	Images
- Wrist Drop - Hand/wrist/finger non-fixed Contracture - MP, PIP, DIP Extension	1. HT-101 **- Hand Thumb 2. RH-101 *- Rest Hand 3. SGH-101 **- Spring Loaded Goniometer Hand	
Post Stroke/CVA (Flaccid/Minimal Tone)	1. HT-101 **- Hand Thumb 2. RH-101 *- Rest Hand 3. OPH-101 *- Opposition Hand Thumb	
Post Stroke/CVA (Spastic/High Tone)	1. F-101 **- Finger Extender 2. C-Grip *- Comfy Grip Hand 3. HAL-101-FS **- Air Hand Roll 4. ACH-101 **- Adjustable Cone Hand	
- Rheumatoid Arthritis - Sprains/Strains/ Fractures/Sports Injury	1. H-101 **- Hand Wrist Finger 2. HT-101 **- Hand Thumb 3. RH-101 *- Rest Hand	
Finger Contractures	1. C-Grip *- Comfy Grip Hand 2. F-101 **- Finger Extender 3. ACH-101 **- Adjustable Cone Hand	
Carpal Tunnel Syndrome	1. Dorsh-101 *- Dorsal Hand 2. F-101(w/o Rolls) **- Finger Extender 3. C-Grip *- Comfy Grip Hand	
Ulnar or Radial Deviation	1. DHT-101 **- Deviation Hand Thumb 2. DRH-101 *- Deviation Rest Hand 3. DF-101 **- Deviation Finger Extender	
Neuromuscular Joint Weakness	1. HT-101 **- Hand Thumb 2. RH-101 *- Rest Hand 3. SGH-101 **- Spring Loaded Goniometer Hand 4. ACH-101 **- Adjustable Cone Hand	
Decreased ROM	1. SGH-101 **- Spring Loaded Goniometer Hand 2. ACH-101 **- Adjustable Cone Hand	
- Individual Finger Positioning (Positioning of MP, PIP & DIP Joints) - Dupuytren's contracture - Boutonniere deformity	1. HSF-101-CP *- Hand Separate Finger (Comfyprene)	

		Sample Indication	Recommended Orthosis	Images		
Elbows		Flaccid/Minimal Tone • Post Stroke/CVA • Post-Surgery • Muscle Weakness • Arthritits	E-101 - Elbow	E-101 GE-101		
		Spastic/High Tone • Non-fixed knee contracture • Post Stroke/CVA • Post-Surgery • Muscle Weakness	GE-101 - Goniometer Elbow	SGE-101		
			SGE-101 - Spring Loaded Goniometer Elbow		LE-101	
			LE-101 - Locking Elbow			EH-101
		Hand & Elbow - 1splint 2 joints • Post Stroke/CVA • Post-Surgery • Muscle Weakness	EH-101 - Elbow Hand Orthosis			
Decreased ROM	SGE-101 - Spring Loaded Goniometer Elbow	GE-101 - Goniometer Elbow Orthosis				
		Sample Indication	Recommended Orthosis	Images		
Knees		Flaccid/Minimal Tone • Post Stroke/CVA • Post-Surgery • Muscle Weakness • Arthritits	K-101 - Knee	K-101 GK-101		
		Spastic/High Tone • Non-fixed knee contracture • Post Stroke/CVA • Post-Surgery • Muscle Weakness	GK-101 - Goniometer Knee	SGK-101		
			SGK-101 - Spring Loaded Goniometer Knee		LK-101	
			LK-101 - Locking Knee			
		Decreased ROM	SGK-101 - Spring Loaded Goniometer Knee	GK-101 - Goniometer Knee		
		Sample Indication	Recommended Orthosis	Images		
Ankle / Foot Universal Sizing		Bedbound Status • Non-fixed ankle contracture	C-Boot - Comfy Boot	C-Boot C-Strap		
		Spasticity • Foot drop (plantar flexion) • Heel Pressure Areas • Plantar Fasciitis • Inversion/Eversion	C-Strap - Comfy Boot with Side Straps (straps control plantar flexion)	A-Boot A-Strap		
			SAF-101 - Spring Loaded Ankle Foot		SAF	
		Ambulating Status • Non-fixed ankle contracture • Foot Drop (plantar flexion) • Heel Pressure Areas • Plantar Fasciitis • Inversion/Eversion	A-Boot - Ambulating Boot			
			A-Strap - Ambulating Boot with Side Strap (straps control plantar flexion)			
			SAF-101 - Spring Loaded Ankle Foot			
		Decreased ROM	SAF-101 - Spring Loaded Ankle Foot Dynamic splint to increase ankle ROM			
		Minimal Aplied Pressure • Stress fractures • Ligament and tendon repair	A-Boot - Ambulating Boot			
	A-Strap - Ambulating Boot with Side Straps					